

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-RNF8-023

Anti RNF8

Product type	Primary antibodies
Immunogen	RNF8 (NP_003949.1, 1 a.a. ~ 485 a.a) full-length human protein.
Form	Liquid
	Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit

Immunogen Sequence

MGEPGFFVTGDRAGGRSWCLRRVGMSAGWLLLEDGCEVTVGRGFGVTYQLVSKICPLMIS
RNHCVLKQNPPEGQWTIMDNKSLNGVWLNRRARLEPLRVYSIHQGDYIQLGVPLENKENAEY
EYEVTEEDWETIYPCLSPKNDQMIEKNKELRTKRKFSLDELAGPGAEGPSNLKSKINKVS
CESGQPVKSQGKGEVASTPSDNLDPKLTALEPSKTTGAPIYPGFPKVTEVHHEQKASNSS
ASQRSLQMFKVTMSRILRLKIQMKEKHEAVMNVKKQTQKGNSKKVVQMEQELQDLQSQLC
AEQAQQQARVEQLEKTFQEEEQHLQGLEIAQGEKDLKQQLAQALQEHWALMEELNRSKKD
FEAIQAKNKELEQTKEEKEKMQAQKEEVLSHMNDVLENELQCIICSEYFIEAVTLNCAH
SFCSYCINEWMKRKIECPICRKDIKSKTYSLVLDNCINKMVNNLSSEVKERRIVLIRERK
AKRLF


Gene Alias DKFZp762K1914, OLD35, PNPASE, old-35

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~55.50 kDa using RNF8 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of RNF8 expression in transfected 293T cell line by RNF8 rabbit polyclonal antibody.

Lane 1: RNF8 transfected lysate(55.50 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619