

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-RPAP3-047

Anti RPAP3

Product type	Primary antibodies
Immunogen	RPAP3 (AAH56415.1, 1 a.a. ~ 631 a.a) full-length human protein.
Form	Liquid
	Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit


Immunogen Sequence

MTSANKAIELQLQVKQNAEELQDFMRDLENWEKDIKQKDMELRRQNGVPEENLPPIRNGN
FRKKKKKGKAKESSKKTREENTKNRIKSYDYEAWAKLDVDRILDELDKDDSTHESLSQESE
SEEDGIHVDSQKALVLKEKGKNKYFKQGKYDEAIDCYTKGMDADPYNPVLPTNRASAYFRL
KKFAVAESDCNLAVALNRSYTKAYSRRGAARFALQKLEEAKKDYERVLELEPNNFEATNE
LRKISQALASKENSYPKEADIVIKSTEGERKQIEAQQNKQQAISEKDRGNGFFKEGKYER
AIECYTRGIAADGANALLPANRAMAYLKIQKYEEAEKDCTQAILLDGSYSKAFARRGTAR
TFLGKLNEAKQDFETVLLLEPGNKQAVTELSKIKKKPLKKVIIETGNLIQTIDVPDSTT
AAPENNPNINLANVIAATGTTSKKNSSQDDLFPSTPTPRAKVLKIEEVSDTSSLQPQASL
KQDVCQSYSEKMPIEIEQKPAQFATTVLPPIPANSFQLESDFRQLKSSPDMLYQYLKQIE
PSLYPKLFQKNLDPDVFNQIVKILHDFYIEKEKPLLFIFEILQRLSELKRFDMAVMFMSET
EKKIARALFNHIDKSGLDSSVEELKKRYGG

Gene Alias FLJ21908, spag

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~71.90 kDa using RPAP3 transfected 293T lysate.

Western Blotting (Transfected lysate)




Western Blot analysis of RPAP3 expression in transfected 293T cell line by RPAP3 rabbit polyclonal antibody.

Lane 1: RPAP3 transfected lysate(71.90 KDa).

Lane 2: Non-transfected lysate.

Antibody dilution: 1:1000. Western Blot analysis of RPAP3 expression in Jurkat (Predicted Mw: 71.8660 kDa).

Western Blotting (Cell lysate)




RPAP3 rabbit polyclonal antibody. Western Blot analysis of RPAP3 expression in Jurkat.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of RPAP3 transfected lysate.

Immunoprecipitation



Immunoprecipitation of RPAP3 transfected lysate using RPAP3 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with RPAP3 purified mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619