

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-RPN2-033

Anti RPN2

Product type	Primary antibodies
Immunogen	RPN2 (NP_002942.2, 1 a.a. ~ 631 a.a) full-length human protein.
Form	Liquid
	Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit

**Immunogen
Sequence**

MAPPGSSTVFLLALTIIASTWALTPHYLTKHDVERLKASLDRPFTNLESAFYISIVGLSS
LGAQVPDAKKACTYIRSNLDPNSVDSLFYAAQASQALSGCEISISNETKDLLLLAAVSEDS
SVTQIYHAVAALSGFGLPLASQEALSALTARLSKEETVLATVQALQTASHLSQQADLRSI
VEEIEDLVARLDELGGVYLQFEEGLETTALFVAATYKLMHDVHGTEPSIKEDQVIQLMNAI
FSKKNFESLSEAFSVASAAVLSHNRYHVPVVVPEGSASDTHEQAILRLQVTNVLSQPL
TQATVKLEHAKSVASRATVLQKTSFTPVGDFELNFMNVKFSSGYDFLVEVEGDNRYIA
NTVELRVKISTEVGITNVDLSTVDKDQSIAPKTTRVTYPAKAKGTFIADSHQNFAFFQL
VDVNTGAELTPHQTFVRLHNQKTGQEVVFVAEPDNKNVYKFELDTSERKIEFDSASGTYT
LYLIIGDATLKNPILWNVADVVIKFPEEEAPSTVLSQNLFTPKQEIQLHFREPEKRPPTV
VSNTFTALILSPLLLLFALWIRIGANVSNTFAPSTIIFHLGHAAMLGLMYVYWTQLNMF
QTLKYLAILGSVTFLAGNRMLAQQAVKRTAH


Gene Alias RIBIIR, RPN-II, RPNII, SWP1

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~69.30 kDa using RPN2 transfected 293T lysate

Western Blotting (Transfected lysate)




Western Blot analysis of RPN2 expression in transfected 293T cell line by RPN2 rabbit polyclonal antibody.

Lane 1: RPN2 transfected lysate(69.30 kDa).

Lane 2: Non-transfected lysate.

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~69.30 kDa using RPN2 transfected 293T lysate

Western Blotting (Tissue lysate)



RPN2 rabbit polyclonal antibody. Western Blot analysis of RPN2 expression in human pancreas.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619