

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-RRS1-393

Anti RRS1

Product type	Primary antibodies
Immunogen	RRS1 (NP_055984.1, 1 a.a. ~ 365 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit

Immunogen Sequence

MEGQSVEELLAKAEQDEAEKLQRITVHKELELQFDLGNLLASDRNPPTGLRCAGPTPEAE
LQALARDNTQLLINQLWQLPTERVEEAIVARLPEPTTRLPREKPLPRPRPLTRWQQFARL
KGIRPKKKTNLVWDEVSGQWRRRWGYQRARDDTKEWLIEVPGNADPLEDQFAKRIQAKKE
RVAKNELNRLNRLARAHKMLPSAAGLHPTGHQSKEELGRAMQVAKVSTASVGRFQERLP
KEKVPRGSGKKRKQPLFGDFAAEKKNQLELLRVMNSKKPQLDVTRATNKQMREEDQEEA
AKRRKMSQKGKRKGGRQGPGGKRKGPPSQGGKRKGGLGGKMNSGPPGLGGKRKGGRPG
GKRRK


KIAA0112

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~41.20 kDa using RRS1 transfected 293T lysate

Western Blotting (Transfected lysate)



Western Blot analysis of RRS1 expression in transfected 293T cell line by RRS1 rabbit polyclonal antibody.


Lane 1: RRS1 transfected lysate(41.20 KDa).

Lane 2: Non-transfected lysate.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of RRS1 transfected lysate.

Immunoprecipitation



Immunoprecipitation of RRS1 transfected lysate using RRS1 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with RRS1 purified mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619