



POLYCLONAL ANTIBODY


For research use only. Not for clinical diagnosis.

Catalog No. CNP-SDPR-167

Anti SDPR

Product type	Primary antibodies
Immunogen	SDPR (NP_004648.1, 1 a.a. ~ 425 a.a) full-length human protein..
Form	Liquid
	Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MGEDAAQAEKFQHPGSMDMRQEKPSPPSPMPSSTPSPSLNLGNTEEAIRDNSQVNAVTVLT LLDKLVNMLDAVQENQHKMEQRQISLEGSVKGIQNDLTKLSKYQASTSNTVSKLLEKSRK VSAHTRAVKERMDRQCAQVKRLNHAQLLRNHFVKVLIFQEENEIPASVFKQPVSGAV EGKEELPDENKSLEETLHTVDLSSDDDLPHDEEALEDSAEKVEESRAEKIKRSSLKKVD SLKKAFSRQNIKKMNKLGTKIVSVERREKIKKSLTSNHQKISSGKSSPFKVSPLTFGRK KVREGESHAENETKSEDLPSSSEQMPNDQEEESFAEGHSEASLASALVEGEIAEEAAEKAT SRGSNSGMDSNIDLTIVEDEEEESVALEQAQKVRYEGSYALTSEEAEERSDGDGPVQPAVLQ VHQTS
Gene alias	PS-p68, SDR
Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~47.20 kDa using SDPR transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of SDPR expression in transfected 293T cell line by SDPR rabbit polyclonal antibody.

Lane 1: SDPR transfected lysate(47.20 KDa).
Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619