

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.


Catalog No. CNP-STX11-014

Anti STX-11

Product type	Primary antibodies
Immunogen	STX11 (NP_003755.2, 1 a.a. ~ 287 a.a) full-length human protein.
Form	Liquid
	Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MKDRLAELLDSLKQYDQQFPDGDDEFDSPHEDIVFETDHILESLYRDIRDIQDENQLLVA DVKRLGKQNARFLTSMRRLSSIKRDTNSIAKAIKARGEVIHCKLRAMKELSEAAEAQHGP HSAVARISRAQYNALTLTFQRAMHDYNQAEMKQRDNCKIRIQRQLEIMGKEVSGDQIEDM FEQGWKWDVFSENLLADVKGARAALNEIESRHRRELLRLESRIKRVHFLQMAVLVEKQAD TLNVIELNVQKTVDYTGQAKAQVRKAVQYEEKNPCRTLCCFCCPCLK

Gene Alias	FHL4, HLH4, HPLH4
Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~33.20 kDa using STX11 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of STX11 expression in transfected 293T cell line by STX11 rabbit polyclonal antibody.


Lane 1: STX11 transfected lysate(33.20 KDa).

Lane 2: Non-transfected lysate.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of STX11 transfected lysate.

Immunoprecipitation



Immunoprecipitation of STX11 transfected lysate using STX11 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with STX11 purified mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619