

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.


Catalog No. CNP-STX12-010

Anti STX-12

Product type	Primary antibodies
Immunogen	STX12 (NP_803173.1, 1 a.a. ~ 276 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MSYGPLDMYRNPGPSGPQLRDFSSIIQTCSGNIQRISQATAQIKNLMSQLGTKQDSSKLQ ENLQQQLQHSTNQLAKETNELLKELGSLPLPLSTSEQRQQRLQKERLMNDFSAALNNFQAV QRRVSEKEKESARARAGSRLSAEERQREEQLVSFDSHEEWNQMMSQEQEDEVATEQDLEL IKERETAIRQLEADILDVNQIFKDLAMMIHDQGDLDISIEANVESSEVHVVERATEQLQRA AYYQKKSRRKMCILVLVLSVILILGLIIWLVYKTK
Gene Alias	MGC51957, STX13, STX14

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~31.60 kDa using STX12 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of STX12 expression in transfected 293T cell line by STX12 rabbit polyclonal antibody.


Lane 1: STX12 transfected lysate(31.60 kDa).

Lane 2: Non-transfected lysate.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of STX12 transfected lysate.

Immunoprecipitation



Immunoprecipitation of STX12 transfected lysate using STX12 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with STX12 mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619