

**POLYCLONAL ANTIBODY**

*For research use only. Not for clinical diagnosis.*

**Catalog No. CNP-TIRAP-119**

# Anti TIRAP

**Product type** Primary antibodies

**Immunogen** TIRAP (AAH32474.1, 1 a.a. ~ 257 a.a) full-length human protein.

**Form** Liquid  
Protein A purified

**Volume** 100ul

**Specificity** Human


**Host** Rabbit

**Immunogen Sequence** MASSTSLPAPGSRPKKPLGKMADWFRQTLLKKPKKRPNSPESTSSDASQPTSQDSPLPPS  
LSSVTSPSLPPTHASDSGSSRWKDYDVCVCHSEEDLVAAQDLVSYLEGSTASLRCFLQL  
RDATPGGAIVSELCQALSSSHCRVLLITPGFLQDPWCKYQMLQALTEAPGAEGCTIPLLS  
GLSRAAYPPELRFMYVVDGRGPDGGFRQVKEAVMRYLQTLSWHLLYHGTPEIGVKLETEN  
PCRASDSHKCDKRYRE\*

**Gene Alias** FLJ42305, Mal, wyatt

**Application notes** **WB**  
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~28.27 kDa using TIRAP transfected 293T lysate.

**Western Blotting (Transfected lysate)**



Western Blot analysis of TIRAP expression in transfected 293T cell line by TIRAP rabbit polyclonal antibody.


Lane 1: TIRAP transfected lysate(28.27 KDa).

Lane 2: Non-transfected lysate.

## IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of TIRAP transfected lysate.

### Immunoprecipitation



Immunoprecipitation of TIRAP transfected lysate using TIRAP rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with TIRAP mouse polyclonal antibody.

## Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

*For research use only. Not for clinical diagnosis.*



**COSMO BIO Co., LTD.**

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp)

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619