



POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-USP18-334

Anti USP18


Product type	Primary antibodies
Immunogen	USP18 (CAK54503.1, 1 a.a. ~ 372 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit

Immunogen Sequence
MSKAFGLLRQICQSILAESSQSPADLEEKKEEDSNMKREQPRERPRAWDYPHGLVGLHNI
GQTCCLNSLIQVFVMNVDFTRILKRITVPRGADEQRRSVPFQMLLLEKMQDSRQKAVRP
LELAYCLQKCNVPLFVQHDAAGLYLKLWNLIKDQITDVHLVERLQALYTIRVKDSLICVD
CAMESSRNSSMLTLPLSLFDVDSKPLKTLEDALHCFQPRELSSKSKCFCECGKKTRGK
QVLKLTHTLPQTLTIHLMRFSIRNSQTRKICHSLYFPQSLDFSQILPMKRESCDAEEQSGG
QYELFAVIAHVGMADSGHYCVYIRNAVDGKWFCFNDSNICLVSWEIQCTYGNPNYHWQE
TAYLLVYMKMEC

iSG43, UBP43

Application notes
WB
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~40.92 kDa using USP18 transfected 293T lysate.

Western Blotting (Transfected lysate)




Western Blot analysis of USP18 expression in transfected 293T cell line by USP18 rabbit polyclonal antibody.

Lane 1: USP18 transfected lysate(40.92 KDa).
Lane 2: Non-transfected lysate.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of USP18 transfected lysate.

Immunoprecipitation



Immunoprecipitation of USP18 transfected lysate using USP18 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with USP18 rabbit polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619