

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-VPS39-185

Anti VPS39


Product type	Primary antibodies
Immunogen	VPS39 (NP_056104.2, 1 a.a. ~ 875 a.a) full-length human protein.
Form	Liquid
	Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit

Immunogen Sequence	MHDAFEPVPILEKLPLQIDCLAAWEEWLLVGTKQGHELLLYRIRKDVGCNRFVETLEKSNK NFSKKIQQIHVVSQFKILVSLLENNIYVHDLTLFQQITTVSKAKGASLFTCDLQHTETGE EVLRCMCVAVKKKLQLYFWKDFHEHQGDFSVDPVPSMAWCENSICVGFKRYYLIRVD GKGSIKELFPTGKQLEPLVAPLADGKVAVGQDDLTVVLNEEGICTQKCALNWTDPVAME HQPPIIYAVLPRYVEIRTFEPRLVQSIELQRPRFITSGGSNIIYVASNHFWRLIPVPM ATQIQQLQDKQFELALQLAEMKDDSDSEKQQQIHHIKNLYAFNLFCQKRFDSEMQVFAK LGTDPHTVMGLYPDLLPTDYRKQLQYPNPLPVLSGAELEKAHLALIDYLTQKRSQVLVKKL NDSHQSSSTPLMEGTPTIKSKKKLLQIIDTLLKCYLHTNVALVAPLLRENNHCHIEE SEHVLKKAHKYSELIIYEKKGLHEKALQVLVDQSKKANSPLKGHERTVQYLQHLGTENL HLIFSYSVWVLRDFPEDGLKIFTEDLPEVESLPRDRVLGFLIENFKGLAIPYLEHIIHVV EETGSRFHNCLIQLYCEKVQGLMKEYLLSFPAGKTPVPAGEEEGELGEYRQKLLMFLEIS SYDPPGRLICDFPFDGLLEERALLGRMGKHEQALFIYVHILKDTRMAEEYCHKHYDRNK DGNKDVYLSLLRMYLSPPSIHCLGPIKLELLEPKANLQAALQVLELHHSKLDTTKALNLL PANTQINDIRIFLEKVLLEENAQKKRFNQVLKNNLHAEFLRVQEERILHQQVKCIITEEKV CMVCKKKIGNSAFARYPNGVHVHYFCSKEVNPADT
---------------------------	--

Gene Alias	FLJ21681, FLJ46546, KIAA0770, TLP, VAM6, hVam6p
-------------------	---

Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~100.80 kDa using VPS39 transfected 293T lysate.
--------------------------	--

Western Blotting (Transfected lysate)



Western Blot analysis of VPS39 expression in transfected 293T cell line by VPS39 rabbit polyclonal antibody.


Lane 1: VPS39 transfected lysate(100.80 KDa).

Lane 2: Non-transfected lysate.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of VPS39 transfected lysate.

Immunoprecipitation



Immunoprecipitation of VPS39 transfected lysate using VPS39 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with VPS39 purified mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619