

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-043

Anti SMARCA5

Product type	Primary antibodies
Immunogen	SMARCA5 (AAH23144.1, 1 a.a. ~ 1052 a.a) full-length human protein.
Form	Liquid
	Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit

**Immunogen
Sequence**

MSSAAEPPPPPPESAPSKPAASIASGGSNSSNKGPEGVAAQAVASAASAGPADAEMEE
IFDDASPGKQKEIQEPDPTYEEKMQTDRANRFEYLLKQTELFHFHQPAAQKTPTSPLKM
KPGRPRIKKDEKQNLLSVGDYRHRRTQEEDDELLTESSKATNVCTRFEDSPSYVKWGKL
RDYQVRGLNWLISLYENGINGILADEMGLGKTLQTISLLGYMKHYRNIPGPHMVLVPKST
LHNWMSEFKRWVPTLRSVCLIGDKEQRAAFVRDVLPGEWDCVTSYEMLIKEKSVFKKF
NWRYLVIDEAHRIKNEKSKLSEIVREFKTTNRLLLTGTPNQNNLHELWLLNLLPDVFN
SADDFDSWFDTNCLGDQKLVERLHMVLRPFLRRIKADVEKSLPPKKEVKIYVGLSKMQ
REWYTRILMKDIDILNSAGKMDKMRLNLMQLRKCCNHPYLFDAEPGPPYTDMHLVT
NSGKMVVLDKLLPKLKEQGSRLIFSQMTRVLDILEDYCMWRNYEYCRLDGGQTPHDERQD
SINAYNEPNSTKFVFMSTRAGGLGINLATADVILYSDWNPQVDLQAMDRAHRIGQTK
TVRVRFITDNTVEERIVERAEMKLRLDSIVIQQGRLVDQNLNKIGKDEMLQMRHGATH
VFASKESITDEDIDILGERGAKKTAEMNEKLSKMGESSLRNFTMDTESSVYNFEGEDYR
EKQKIAFTIEWIEPPKRERKANYAVDAYFREARLVSEPKAPKAPRPPKQPNVQDFQFFPPR
LFELLEKEILFYRKTYGYKVPNPPELNAQAQKEEQLKIDEAESLNDEELEEKEKLLTQ
GFTNWNKRDFNQFIKANKEWGRDDIENIAREVEGKTPEEVIEYSVFWERCNELQDIEKI
MAQIERGEARIQRRISIKKALDTKIGRYKAPFHQLRISYGTNKGKNYTEEEDRFLICMLH
KLGFDKENVYDELRCIRNSPQFRFDWFLKSRTAMELQRRCNTLITLIERENMELEEKEK
AEKKKRGPKPSTQKRKMDGAPDGRGRKKKLKL


Gene Alias AMAC-1, AMAC1, CKb7, DC-CK1, DCCCK1, MIP-4, PARC, SCYA18

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~121.90 kDa using SMARCA5 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of SMARCA5 expression in transfected 293T cell line by SMARCA5 rabbit polyclonal antibody.

Lane 1: SMARCA5 transfected lysate(121.90 kDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619