



POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-058

Anti CCL14


Product type	Primary antibodies
Immunogen	CCL14 (NP_004157.1, 1 a.a. ~ 93 a.a) full-length human protein.
Form	Liquid
	Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit

Immunogen Sequence	MKISVAaipffllitlalgktesssrgpyhpsccftyttykiprqrimdyyetnsqcs kpgivfitkrghsvctnpsdkwvqdyikdmken
---------------------------	---

Gene Alias	CC-1, CC-3, CKb1, HCC-1, HCC-3, MCIF, NCC-2, NCC2, SCYA14, S CYL2, SY14
-------------------	--

Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~10.70 kDa using CCL14 transfected 293T lysate.
--------------------------	---

Western Blotting (Transfected lysate)



Western Blot analysis of CCL14 expression in transfected 293T cell line by CCL14 rabbit polyclonal antibody.

Lane 1: CCL14 transfected lysate(10.70 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619