

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-068


Anti UBE3C

Product type	Primary antibodies
Immunogen	MAPKAPK2 (NP_116584.2, 1 a.a. ~ 400 a.a) full-length human protein.
Form	Liquid
	Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit

Immunogen Sequence	MLSNSQGQSPPVFPAPAPPPQPPTPALPHPPAQPPPPPPQFPQFHVKSGLQIKKNAII DDYKVTSQVLGLGINGKVLQIFNKRTQEFALKMLQDCPKARREVELHWRASQCPIHVR VDVYENLYAGRKCLLIVMECLDGGELFSRIQDRGDQAFTEREASEIMKSIGEGAIQYLHSI NIAHRDVKPENLLYTSKRPNAILKLTDFGFAKETTSNLSLTTPCYTPYYVAPEVLGPEKY DKSCDMWSLGVIMYILLCGYPPFYSNHGLAISPGMKTRIRMGQYEFNPPEWSEVSEEVKM LIRNLLKTEPTQRMTEFMNHPWIMQSTKVPQTPLHTSRVLKEDKERWEDVKEEMTSAL ATMRVDYEQIKIKKIEDASNPLLLKRRKKARALEAAALAH
---------------------------	--

Gene Alias	UBE3C (AAH26241.1, 1 a.a. ~ 404 a.a) full-length human protein.
Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~45.80 kDa using UBE3C transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of UBE3C expression in transfected 293T cell line by UBE3C rabbit polyclonal antibody.

Lane 1: UBE3C transfected lysate(45.80 KDa).
Lane 2: Non-transfected lysate.

I

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619