

**POLYCLONAL ANTIBODY**

*For research use only. Not for clinical diagnosis.*

**Catalog No. CNP- WB-094**

# Anti ADHFE1

**Product type** Primary antibodies

**Immunogen** ADHFE1 (NP\_653251.1, 1 a.a. ~ 419 a.a) full-length human protein.

**Form** Liquid

Protein A purified

**Volume** 100ul

**Specificity** Human

**Host** Rabbit


**Immunogen Sequence** MAVSNIRYGAAVTKEVGMDLKNMGAKNVCLMTDKNLSKLPPVQVAMDSLVKNGIPFTVYD  
NVRVEPTDSSFMEAEIFAQKGAFDAYVAVGGGSTMDTCKAANLYASSPHSDFLDYVSAPI  
GKGKPVSVPLKPLIAVPPTTSGTGSETTGVAIFDYEHLKVKIGITSRAIKPTLGLIDPLHT  
LHMPARVVANSFGFDVLCHALESYTTLPYHLRSPCPSNPITRPAYQGSNPISDIWAIHALR  
IVAKYLKRAVRNPDDLEARSHMHLASAFAGIGFGNAGVHLCHGMSYPISGLVKMYKAKDY  
NVDHPLVPHGLSVVLTSPAVFTFTAQMFPERHLEMAEILGADTRTARIQDAGLVLADTLR  
KFLFDLDVDDGLAAGVGYSKADIPALVKGTLQPQERVTKLAPRPQSEEDLAALFEASMKLY

**Gene Alias** ADH8, FLJ32430, HMFT2263, HOT, MGC48605

**Application** **WB**

**notes** Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~45.20 kDa using ADHFE1 transfected 293T lysate.

**Western Blotting (Transfected lysate)**




Western Blot analysis of ADHFE1 expression in transfected 293T cell line by ADHFE1 rabbit polyclonal antibody.

Lane 1: ADHFE1 transfected lysate(45.20 KDa).

Lane 2: Non-transfected lysate.

Antibody dilution: 1:1000. Western Blot analysis of ADHFE1 expression in K-562 (Predicted Mw: 50.3080 kDa).


**Western Blotting (Cell lysate)**



ADHFE1 rabbit polyclonal antibody. Western Blot analysis of ADHFE1 expression in K-562.

Antibody dilution: 1:1000. Western Blot analysis of ADHFE1 expression in mouse kidney (Predicted Mw: 50.3080 kDa).

**Western Blotting (Tissue lysate)**



ADHFE1 rabbit polyclonal antibody. Western Blot analysis of ADHFE1 expression in mouse kidney.

**Storage**

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

*For research use only. Not for clinical diagnosis.*



**COSMO BIO Co., LTD.**

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp)

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619