

POLYCLONAL ANTIBODY


For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-118

Anti ZNF92

Product type	Primary antibodies
Immunogen	ZNF92 (AAH36439.1, 1 a.a. ~ 517 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MVDKTPVMCSHFAQDVWPEHSIKDSFQKVILRTYGKYGHENLQLRKDHKSVDACKVYKGG YNGLNQCLTTTDSKIFQCDKYVKVFHKFPNVNRNKIRHTGKKPFKCKNRGKSFCMLSQLT QHKKIHTREYSYKCEECGKAFNWSSTLTCHKIHTGEKPYKCEECGKAFNRSSNLTKHKI IHTGEKPYKCEECGKAFNRSSLTCHKRIHTEEKPYKCEECGKAFNQFSILNKHKRIHME DKPYKCEECGKAFRVFSILKHKHIIHTGEKPYKCEECGKAFNQFSNLTKHKIHTGEKPY KCDECGKAFNQSSLTCHKRIHTGEKPYKCEECGKAFKQSSTLTEHKIHTGEKPYKCEK CGKAFSWSSAFTKHKRNHMDKPYKCEECGKAFSVFSTLTCHKIHTREKPYKCEECGKA FNRSSIFTKHKIHTTEGKSYKCEKCGNAFNQSSNLTARKIIYTGEKPYKYEECDKAFNKF STLITHQIIYTGEKPCKHECGRAFKNSSNYTKEKLQT
Gene Alias	HPF12, HTF12, TF12
Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~60.20 kDa using ZNF92 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of ZNF92 expression in transfected 293T cell line by ZNF92 rabbit polyclonal antibody.

Lane 1: ZNF92 transfected lysate(60.20 kDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619