

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-136

Anti IGSF4B

Product type Primary antibodies

Immunogen IGSF4B (AAH33819.1, 1 a.a. ~ 398 a.a) full-length human protein

Form Liquid

Protein A purified

Volume 100ul

Specificity Human


Host Rabbit

Immunogen Sequence MGAPAAASLLLLLLFACCWAPGGANLSQDDSQPWTSDET VVAGGT VVLKCQVKDHEDSSL
QWSNPAQQTL YFGEKRALRDNRIQLVTSTPHEL SISISNVALADEGEYTCSIFTMPVRTA
KSLVTVLGIPQKPIITGYKSSLREKDTATLNCQSSGSKPAARLTWRKGDQELHGEPTRIQ
EDPNGKTFTVSSSVTFQVTREDDGASIVCSVNHESLKGADRSTSQRIEVLYTPTAMIRPD
PPHPREGQKLLHCEGRGNPVPQQYLWEKEGSPPLKMTQESALIFPFLNKSDSGTYGCT
ATSNMGSYKAYYTLNVNDPSPVPSSSSTYHAIIGGIVAFIVFLLIMLIFLGHYLIRHKG
TYLTAEAKGSDDAPDADTAIINAEGGQSGGDDKKEYFI

Gene Alias GCP2, MGC138162, SPBC97, Spc97p

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~43.30 kDa using CADM3 transfected 293T lysate..

Western Blotting (Transfected lysate)



Western Blot analysis of CADM3 expression in transfected 293T cell line by CADM3 rabbit polyclonal antibody

Lane 1: CADM3 transfected lysate(43.30 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp> e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617 [国内連絡先] Phone : +81-3-5632-9610
FAX : +81-3-5632-9618 F AX : +81-3-5632-9619