

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-143

Anti P2RX1

Product type Primary antibodies

Immunogen P2RX1 (NP_002549.1, 1 a.a. ~ 399 a.a) full-length human protein.

Form Liquid

Protein A purified

Volume 100ul

Specificity Human

Host Rabbit

Immunogen Sequence MARRFQEELAAFLFEYDTPRMVLVRNKKVGVIFRLIQLVVLVYVIGWVFLYEKGYQTSSG
LISSVSVKLKGLAVTQLPGLGPQVWDVADYVFPAQGDNSFVVMTNFIVTPKQTQGYCAEH
PEGGICKEDSGCTPGKAKRKAQGIRTGKCVAFNDTVKTCEIFGWCPVEVDDDDIPRPALLR
EAENFTLFIKNSISFPRFKVNRNRLVEEVNAAHMKTCLFHKTLLHPLCPVFQLGYVYVQESG
QNFSTLAEKGGVVGITIDWHCDLDWHVRHCRPIYEFHGLYEEKNLSPGFNFRFARHFVEN
GTNYRHLFKVFGIRFDILVDGKAGKFDIIPMTTIGSGIGIFGVATVLCDLLLHILPKR
HYYKQKKFKYAEDMGPGAAERDLAATSSTLGLQENMRTS


Gene Alias P2X1

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~45.00 kDa using P2RX1 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of P2RX1 expression in transfected 293T cell line by P2RX1 rabbit polyclonal antibody.

Lane 1: P2RX1 transfected lysate(45.00 kDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619