



POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.


Catalog No. CNP- WB-163

Anti CLEC9A

Product type	Primary antibodies
Immunogen	CLEC9A (NP_997228.1, 1 a.a. ~ 241 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MHEEEIYTSLQWDSPAPDTYQKCLSSNKCSGACCLVMVISCVFCMGLLTASIFLGVKLLQ VSTIAMQQQEKLIIQKERALLNFTIEWKRSCALQMKYCQAFMQNSLSSAHNSSPCPNNWIQN RESCYYVSEIWSIWHTSQENCLKEGSTLLQIESKEEMDFITGSLRKIKGSYDYWVGLSQD GHSGRWLWQDGSSPSGLLPAERSQSANQVCGYVKSNSLLSSNCSTWKYFICEKYALRSS V

Gene Alias	DNGR1, UNQ9341
Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~27.30 kDa using CLEC9A transfected 293T lysate.

Western Blotting (Transfected lysate)




Western Blot analysis of CLEC9A expression in transfected 293T cell line by CLEC9A rabbit polyclonal antibody.

Lane 1: CLEC9A transfected lysate(27.30 KDa).

Lane 2: Non-transfected lysate.

Antibody dilution: 1:1000. Western Blot analysis of CLEC9A expression in human spleen (Predicted Mw: 27.3240 kDa).

Western Blotting (Tissue lysate)



CLEC9A rabbit polyclonal antibody. Western Blot analysis of CLEC9A expression in human spleen.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619