



POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-183

Anti CCRL2

Product type Primary antibodies

Immunogen CCRL2 (NP_003956.2, 1 a.a. ~ 344 a.a) full-length human protein.

Form Liquid
Protein A purified

Volume 100ul

Specificity Human

Host Rabbit

Immunogen Sequence MANYTLAPEDEYDVLIEGELESDEAEQCDKYDAQALSAQLVPSLCSAVFVIGVLDNLLVV
LILVKYKGLKRVENIYLLNLAVSNLCFLTLFPWAHAGGDPMCKILIGLYFVGLYSETFF
NCLLTVQRYLVFLHKGNNFSARRRVPCGIITSVLAWVTAILATLPEFVVYKPKMEDQKYK
CAFSRTPFLPADETFWKHFLTLKMNISVLVPLFIFTFLYVQMRKTLRFREQRYSLFKLV
FAIMVVFLLMWAPYNIAFFLSTFKEHFSLSDCKSSYNLDKSVHITKLIATTHCCINPLLY
AFLDGTFSKYLRCRCFHLRSNTPLQPRGQSAQGTSREEPDHSTEV


Gene Alias CKRX, CRAM-A, CRAM-B, FLJ55815, HCR, MGC116710

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~39.50 kDa using CCRL2 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of CCRL2 expression in transfected 293T cell line by CCRL2 rabbit polyclonal antibody.

Lane 1: CCRL2 transfected lysate(39.50 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619