

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.


Catalog No. CNP- WB-219

Anti DLL1

Product type	Primary antibodies
Immunogen	DLL1 (NP_005609.3, 1 a.a. ~ 723 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MGSRCALALAVLSALLCQVWSSGVFELKLQEFVNKKGLLGNRNCCRGGAGPPPCACRTFF RVCLKHYQASVSPEPPCTYGSVATPVLGVDSESLPDGGGADSAFSNPIRFPFGFTWPGTF SLIIEALHTDSPDDLATENPERLISRLATQRHLTVGEEWSQDLHSSGRTDLKYSYRFVCD EHYYGEGCSVFCRPRDDAFGHFTCGERGEKVCNPGWKGPYCTEPICLPGCDEQHGFCDKP GECKCRVGWQGRYCDECIRYPGCLHGTCQQPWQCNCQEGWGGLFCNQDLNYCTHHKPCKN GATCTNTGQGSYTCSCRPGYTGATCELGIDECDPSPCKNGGSCDLENSYSCTCPPGFY KICELSAMTCADGPCFNGGRCSDSPDGGYSCRCVGYSGFNCEKKIDYCSSSPCSNGAKC VDLGDAYLCRCQAGFSGRHCDNDVDDCASSPCANGGTCRDGVNDFSTCPPGYTGRNCSA PVSRCEHAPCHNGATCHERGHRYVCECARGYGGPNCQFLLPELPPGPAVVDLTKLEGQG GPFPWVAVCAGVILVLMLLGCAAVVVCVRLRLQKHRPPADPCRGETETMNNLANCQREK DISVSIIGATQIKNTNKKADFGDHSADKNGFKARYPAVDYNLVQDLKGDDTAVRDAHSK RDTKCQPQGSSGEEKGTPTTLRGGEASERKRPDSGCSTSKDTKYQSVYVISEEKDECVIA TEV
Gene Alias	DELTA1, DL1, Delta

Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~79.53 kDa using DLL1 transfected 293T lysate.
--------------------------	---

Western Blotting (Transfected lysate)




Western Blot analysis of DLL1 expression in transfected 293T cell line by DLL1 rabbit polyclonal antibody.

Lane 1: DLL1 transfected lysate(79.53 KDa).

Lane 2: Non-transfected lysate.

Antibody dilution: 1:1000. Western Blot analysis of DLL1 expression in Jurkat (Predicted Mw: 77.9560 kDa).

Western Blotting (Cell lysate)



DLL1 rabbit polyclonal antibody. Western Blot analysis of DLL1 expression in Jurkat.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619