

**POLYCLONAL ANTIBODY**


*For research use only. Not for clinical diagnosis.*

**Catalog No. CNP- WB-220**

# Anti CHKA

|                           |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
|---------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product type</b>       | Primary antibodies                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Immunogen</b>          | CHKA (NP_001268.2, 1 a.a. ~ 457 a.a) full-length human protein.                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Form</b>               | Liquid<br>Protein A purified                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Volume</b>             | 100ul                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Specificity</b>        | Human                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Host</b>               | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Immunogen Sequence</b> | MKTKFCTGGEAEPSPLGLLSCGSGSAAPAGVGQQRDAASDLESKQLGGQQPPLALPPP<br>PPLPLPLPLPQPPPPQPPADEQPEPRTRRRAYLWCKEFLPGAWRGLREDEFHISVIRGGL<br>SNMLFQCSLPDATTATLGDEPRKVLLRLYGAILQMRSCNKEGSEQAQKENEFGAEAMVLE<br>SVMFAILAERSLGPPLYGIFPQGRLEQFIPSRRLDTEELSLPDIAEIAEKMATFHGMKM<br>PFNKEPKWLFGTMEKYLKEVLRKFTEESRIKKLHKLLSYNLPLELENLRSLESTPSPV<br>VFCHNDCQEGNILLLEGRENSEKQKLMLIDFEYSSYNRYRGFDIGNHFCEWMYDYSYEKYP<br>FFRANIRKYPTKKQQLHFISYLPAPFQNDNFENLSTEEKSIIKEEMLLEVNRALASHFLW<br>GLWSIVQAKISSIEFGYMDYAQARFDAYFHQKRKLGV |
| <b>Gene Alias</b>         | CHK, CKI                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Application notes</b>  | <b>WB</b><br>Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~52.20 kDa using CHKA transfected 293T lysate.                                                                                                                                                                                                                                                                                                                                                                                  |

**Western Blotting (Transfected lysate)**



Western Blot analysis of CHKA expression in transfected 293T cell line by CHKA rabbit polyclonal antibody.

Lane 1: CHKA transfected lysate(52.20 KDa).  
Lane 2: Non-transfected lysate.

**Storage**

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

*For research use only. Not for clinical diagnosis.*



**COSMO BIO Co., LTD.**

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp)

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619