

POLYCLONAL ANTIBODY


For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-221

Anti DBT

Product type	Primary antibodies
Immunogen	DBT (AAH16675.1, 1 a.a. ~ 482 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MAAVRMLRTWSRNAGKLCVRYFQTCGNVHVLKPNYVCFFGYPSFKYSHPHHFLKTTAAL RGQVVQFKLSDIGEGIREVTVKWYVKEGDTVVSQFDSICEVQSDKASVTITSRYDGVKK LYYNLDDIAYVGKPLVDIETEALKDSEEDVVETPAVSHDEHHTQEIKGRKTLATPAVRRRL AMENNIKLSSEVVGSGKDGRILKEDILNYLEKQTGAILPPSPKVEIMPPPPKPKDMTVPIIL VSKPPVFTGKDKTEPIKGFQKAMVKTMSAALKIPHFGYCDEIDLTELVLKREELKPIAFA RGIKLSFMPFFLKAASLGLLQFPILNASVDENCQNITYKASHNIGIAMDTEQGLIVPNVK NVQICSIFDIATELNRLQKLGSVGLSTTDLTGGTFTLSNIGSIGGTFAKPVIMPPEVAI GALGSIKAIPRFNQKGEVYKAQIMNVSWWSADHRVIDGATMSRFSNLWKSYLENPAFMLLD LK
Gene Alias	BCATE2, E2, E2B, MGC9061
Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~53.50 kDa using DBT transfected 293T lysate.

Western Blotting (Transfected lysate)




Western Blot analysis of DBT expression in transfected 293T cell line by DBT rabbit polyclonal antibody.

Lane 1: DBT transfected lysate(53.50 KDa).
Lane 2: Non-transfected lysate.

Antibody dilution: 1:1000. Western Blot analysis of DBT expression in A-431 (Predicted Mw: 53.4870 kDa).


Western Blotting (Cell lysate)



DBT rabbit polyclonal antibody. Western Blot analysis of DBT expression in A-431.

Antibody dilution: 1:1000. Western Blot analysis of DBT expression in human kidney (Predicted Mw: 53.4870 kDa).

Western Blotting (Tissue lysate)



DBT rabbit polyclonal antibody. Western Blot analysis of DBT expression in human kidney.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619