

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-223


Anti MANBA

Product type	Primary antibodies
Immunogen	MANBA (AAH15743.1, 1 a.a. ~ 879 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MRLHLLLLLALCGAGTTAAELSYSLRGNWSICNGNGSLELPGAVPGCVHSALFQQGLIQD SYRFDNLNLRWVSLDNWTYSKEFKIPFEISKWQKVNLIIEGVDTVSKILFNEVTIGETD NMFNRYSFIDITNVVRDVNSIELRFQSAVLYAAQQSKAHTRYQVPPDCPPLVQKGECHVNF VRKEQCSFSWDWGPSFPTQGIWKDVRIEAYNICHNLNYFTFSPIYDKSAQEWNLEIESTFD VVSSKPVGGQVIVAIPKLQTQQTYSIELQPGKRIVELFVNISKNITVETWWPHGHGNQTG YNMTVLFELDGGLNIEKSAKVYFRTVELIEEPIKGSPGLSFYFKINGFPIFLKGSNWIPA DSFQDRVTSELLRLLLQSVVDANMNTLRVWGGGIYEQDEFYELCDELGIMVWQDFMFACA LYPTDQGF LDSVTAEVAYQIKRLKSHPSIIWSGNNEEALMMNWHISFTDRPIYIKD YVTLYVKNIRELVLAGDKSRPFITSSPTNGAETVAEAWVSQNPNSNYFGDVHFYDYISDC WNWKVFPKARFASEYGYQSWPSFSTLEKVSSTEDWSFNKFSLHRQHHEGGNKQMLYQAG LHFKLPQSTDPLRTFKDTIYLTQVMQAQC VKTETEFYRRSRSEIVDQQGHTMGALYWQLN DIWQAPSWASLEYGGKWKMLHYFAQNFFAPLLPVGFENENMFYIYGVS DLHSDYSMTLSV RVHTWSSLEPVCSRVTFRFVMKGGEAVCLYEPPVSELLRRCGNCTRESCVVSFYLSADHE LLSPTNYHFLSSPKEAVGLCKAQITAIISQQGDIFVFDLETSAPAPFVWLDVGSIPGRFS DNGFLMTEKTRTILFYPWEPTSKNELEQSFHVTS LTDIY

Gene Alias MANB1

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~100.90 kDa using MANBA transfected 293T lysate

Western Blotting (Transfected lysate)



Western Blot analysis of MANBA expression in transfected 293T cell line by MANBA rabbit polyclonal antibody.

Lane 1: MANBA transfected lysate(100.90 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619