

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-240

Anti POLE

Product type Primary antibodies

Immunogen POLE (AAH21559.1, 1 a.a. ~ 370 a.a) full-length human protein.

Form Liquid
Protein A purified

Volume 100ul

Specificity Human

Host Rabbit


Immunogen Sequence MDPSNYGGIKGKVSSRIHCGLQDSQKAGGAEDEQENEDDEEERDGEAAAAEESNVEDLL
ENNWNILQFLPQAASCQNYFLMIVSAYIVAVYHCMKDGLRRSAPGSTPVRRRGASQLSQE
AEGAVGALPGMITFSQDYVANELTQSFFTITQKIQKKVTGSRNSTELSEMFPVLPGSHLL
LNNPALEFIKYVCKVLSLDTNITNQVNKLNRDLLRLVDVGEFSEEAQFRDPCRSYVLPEV
ICRSCNFCRDLDLCKDSSFSEDGAVLPQWLCSNCQAPYDSSAIEMTLVEVLQKKLMAFTL
QDLVCLKCRGVKETSMPVYCSAGDFALTHTQVFMEQIGIFRNIAQHYGMSYLLETLEW
LLQKNPQLGH

Gene Alias DKFZp434F222, FLJ21434, POLE1

Application **WB**

notes Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~41.50 kDa using POLE transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of POLE expression in transfected 293T cell line by POLE rabbit polyclonal antibody.

Lane 1: POLE transfected lysate(41.50 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619