

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-243

Anti MAN1A1

Product type Primary antibodies

Immunogen MAN1A1 (AAH65827.1, 1 a.a. ~ 367 a.a) full-length human protein.

Form Liquid
Protein A purified

Volume 100ul

Specificity Human


Host Rabbit

Immunogen Sequence MPVGGLPLFSSPAGGVLGGGLGGGGGRKGS GPAALRLTEKFVLLLVFSAFITLCFGAIF
FLPDSSKLLSGVLFHSSPALQPAADHKPGPGARAEDAAEGRARRREEGAPGDPEAALEDN
LARIRENHERALREAKETLQKLPEEIQRDILLEKKKVAQDQLRDKAPFRGLPPVDFVPPI
GVESREPADAAIREKRAKIKEEIEARKGQANCPGSSSSCEVEIQMMKHAWNNYKGYAWGL
NELKPISKGGHSSSLFGNIKGATIVDALDTLFIMEMKHEFEEAKSWVEENLDFNVNAEIS
VFEVNIRFVGGLLSAYYLSGEEIFRKKAVELGVKLLPAFHTPSGIPWALLNMKRWMCNIP
RAMEATI

Gene Alias HUMM3, HUMM9, MAN9

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~40.00 kDa using MAN1A1 transfected 293T lysate.

Western Blotting (Transfected lysate)




Western Blot analysis of MAN1A1 expression in transfected 293T cell line by MAN1A1 rabbit polyclonal antibody.

Lane 1: MAN1A1 transfected lysate(40.00 KDa).

Lane 2: Non-transfected lysate.

Antibody dilution: 1:1000. Western Blot analysis of MAN1A1 expression in A-431 (Predicted Mw: 72.9690 kDa).

Western Blotting (Cell lysate)



MAN1A1 rabbit polyclonal antibody. Western Blot analysis of MAN1A1 expression in A-431.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619