

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-261

Anti SGPL1

Product type Primary antibodies

Immunogen SGPL1 (NP_003892.2, 1 a.a. ~ 568 a.a) full-length human protein.

Form Liquid
Protein A purified

Volume 100ul

Specificity Human


Host Rabbit

Immunogen Sequence MPSTDLLMLKAFEPYLEILEVYSTKAKNYVNGHCTKYEPWQLIAWSVVWTLIVWGYEFV
FQPESLWSRFKKKCFKLTRKMPIIGRKIQDKLNKTKDDISKNMSFLKVDKEYVKALPSQG
LSSSAVLEKLKEYSSMDAFWQEGRASGTVYSGEEKL TELLVKAYGDFAWSNPLHPDIFPG
LRKIEAEIVRIACSLFNGGPDSCGCVTSGGTESILMACKAYRDLAFEKGIKTPEIVAPQS
AHAAFNKAASYFGMKIVRVPLTKMMEVDVVRAMRRAISRNTAMLCSTPQFPHGVIDPVE
VAKLAVKYKIPLHVDACLGGFLIVFMEKAGYPLEHPDFRVKGVTSISADTHKYGYAPKG
SSLVLYSDKKYRNYQFFVDTDWQGGIYASPTIAGSRPGGISAACWAALMHFGENGYVEAT
KQIIKTARFLKSELENIKGIFVFGNPQLSVIALGSRDFDIYRLSNLMTAKGWNLNQLQFP
PSIHFCITLLHARKRVAIQFLKDIRESVTQIMKNPKAKTTGMGAIYGMAQTTVDNRNMVAE
LSSVFLDSLSTDTVTQGSQMNGSPKPH

Gene Alias FLJ13811, KIAA1252, SPL

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~63.50 kDa using SGPL1 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of SGPL1 expression in transfected 293T cell line by SGPL1 rabbit polyclonal antibody.

Lane 1: SGPL1 transfected lysate(63.50 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619