

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-274


Anti PPA2

Product type	Primary antibodies
Immunogen	PPA2 (NP_789845.1, 1 a.a. ~ 334 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MSALLRLLRTGAPAAACLRRLGTSAGTGSRRAMALYHTEERGQPCSQNYRLFFKNVTGHI SPFHDIPLKVNSKEENGIPMKKARNDYENLFNMIVEIPRWTNAKMEIATKEPMNPIKQY VKDGKRLRYVANIFPYKGYIWNYGTLPTQWEDPHEKDKSTNCFGDNDPIDVCEIGSKILSC GEVIHVKILGILALIDEGETDWKLIANANDPEASKFHDIDDVKKFKPGYLEATLNWFRL YKVPDGKPENQFAFNGEFKNKAFALVIVKSTHQCWKALLMKKCNGGAINCTNVQISDSPF RCTQEEARSLVESVSSSPNKESNEEEQVWHFLGK

Gene Alias FLJ20459, HSPC124, MGC49850, SID6-306

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~37.90 kDa using PPA2 transfected 293T lysate.

Western Blotting (Transfected lysate)




Western Blot analysis of PPA2 expression in transfected 293T cell line by PPA2 rabbit polyclonal antibody.

Lane 1: PPA2 transfected lysate(37.90 KDa).

Lane 2: Non-transfected lysate.

Antibody dilution: 1:1000. Western Blot analysis of PPA2 expression in rat brain (Predicted Mw: 37.9620 kDa).

Western Blotting (Tissue lysate)



PPA2 rabbit polyclonal antibody. Western Blot analysis of PPA2 expression in rat brain.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619