



POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.


Catalog No. CNP- WB-275

Anti ACAD8

Product type	Primary antibodies
Immunogen	ACAD8 (AAH01964.1, 1 a.a. ~ 415 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MLWSGCRRFGARLGCLPGGLRVLVQTGHRSLTSCIDPSMGLNEEQKEFQKVAFDFAAREM APNMAEWDQKELFPVDVMRKAQGLFGGVYIQTDVGGSGLSRLDTSVIFEALATGCTSTT AYISIHNMCAWMIDSFSGNEEQRHKFCPLCTMEKFASYCLTEPGSGSDAASLLTSAKKQG DHYILNGSKAFISGAGESDIYVVMCRTGGLGPKGISCIVVEKGTPGLSFGKKEKKVGWNS QPTRAVIFEDCAVPVANRIGSEGQGFLIAVRGLNGGRINIASCSLGAAHASVILTRDHLN VRKQFGEPLASNQYLQFTLADMATRLVAARLMVRNAAVALQEERKDAVALCSMAKLFATD ECFAICNQALQMGGYGYLKDYAVQQYVRDSRVHQILEGSNEVMRILISRLLQE

Gene Alias	ACAD-8, FLJ22590
Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~45.10 kDa using ACAD8 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of ACAD8 expression in transfected 293T cell line by ACAD8 rabbit polyclonal antibody.

Lane 1: ACAD8 transfected lysate(45.10 KDa).
Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619