


POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-287

Anti GNS

Product type	Primary antibodies
Immunogen	GNS (NP_002067.1, 1 a.a. ~ 552 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MRLPLAPGRLRRGSPRHLPS CSPALLLVGGCLGVFGVAAGTRRPNVVLLLTDDQDEV LGGMTPLKKTALIGEMGMTFSSAYVPSALCCPSRASILTGKYPHNHHVNNNTLEGNCSS KSWQKIQEPNTFPAILRSMCGYQTFFAGKYLNEYGAPDAGGLEHVPLGWSYWYALEKNSK YYNYTSLINGKARKHGENYSVDYLTDLANVSLDFLDYKSNFEPFFMMIATPAPHSPWTA APQYQKAFQNVFAPRNKNFNIHGTNKHWRIRQA KTPMTNSSIQFLDNAFRKRWQTLLSVD DLVEKLVKRLEFTGELNNTYIFYTSDNGYHTGQFSLPIDKRQLYEFDIKVPLLVRGPGIK PNQTSKMLVANIDLGPTILDIAGYDLNKTQMDGMSLLPILRGASNLTWRS DVLVEYQGEG RNVTDPTCPSLSPGV SQCFPDCVCEDAYNNTYACVRTMSALWNLQYCEFD DQEVFVEVYN LTADPDQITNIAKTIDPELLGKMNYRLMMLQSCSGPTCRTPGVFD PGYRFDPRLMFSNRG SVRTRRFSKHL
Gene Alias	G6S, MGC21274
Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~62.10 kDa using GNS transfected 293T lysate. Western Blotting (Transfected lysate) 


Western Blot analysis of GNS expression in transfected 293T cell line by GNS rabbit polyclonal antibody.

Lane 1: GNS transfected lysate(62.10 KDa).

Lane 2: Non-transfected lysate.

Antibody dilution: 1:1000. Western Blot analysis of GNS expression in human liver (Predicted Mw: 62.0820 kDa).

Western Blotting (Tissue lysate)



GNS rabbit polyclonal antibody. Western Blot analysis of GNS expression in human liver.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619