

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-296

Anti HS3ST3A1

Product type Primary antibodies

Immunogen HS3ST3A1 (NP_006033.1, 1 a.a. ~ 406 a.a) full-length human protein.

Form Liquid
Protein A purified

Volume 100ul

Specificity Human


Host Rabbit

Immunogen Sequence MAPPGPASALSTSAEPLSRSIFRKFLMLCSLLTSLYVFYCLAERCQTLSGPVVGLSGGG
EEAGAPGGGVLGGPRELAVWPAAAQRKRLQLPQWRRRRPPAPRDDGEEAAWEEESPGL
SGGPGGSGAGSTVAEAPPGLALLLDEGSKQLPQAIIGVKKGGTRALLEFLRVHPDVRA
VGAEPHFDFRSYDKGLAWYRDLMPRTLDDGQITMEKTPSYFVTREAPARISAMSKDTKLIV
VVRDPVTRAISDYTQTLKRDPDIPTFESLTFKNRTAGLIDTSWSAIQIGIYAKHLEHWLR
HFPIRQMLFVSGERLISDPAGELGRVQDFLGLKRIITDKHFYFNKTKGFPCLKKAEGSSR
PHCLGKTKGRTHPEIDREVVRRLREFYRPFNLKFYQMTGHDFGWDG

Gene Alias 3OST3A1, 3OST3A1

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~44.90 kDa using HS3ST3A1 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of HS3ST3A1 expression in transfected 293T cell line by HS3ST3A1 rabbit polyclonal antibody.

Lane 1: HS3ST3A1 transfected lysate(44.90 KDa).
Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619