



POLYCLONAL ANTIBODY


For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-298

Anti POLR3F

Product type	Primary antibodies
Immunogen	POLR3F (NP_006457.2, 1 a.a. ~ 316 a.a) full-length human protein..
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MAEVKVKVQPPDADPVEIENRIIELCHQFPHGITDQVIQNEMPHIEAQQRAVAINRLLSM GQLDLLRSNTGLLYRIKDSQNAGKMKGSDNQEKLVIYQIIEDAGNKGIVSRDIRYKSNLPL TEINKILKNLESKKLIKAVKSVAASKKKVYMLYNLQPDRSVTGGAWYSDQDFESEFVEVL NQQCFKFLQSKAETARESKQNPMIQRNSSFASSHEVWKYICELGISKVELSMEDIETILN TLIYDGKVEMTIIAAKEGTVGSDGHMKLYRAVNPIIPPTGLVRAPCGLCPVFDDCHEGG EISPSNCIYMTEWLEF
Gene Alias	MGC13517, RPC39, RPC6
Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~35.70 kDa using POLR3F transfected 293T lysate.

Western Blotting (Transfected lysate)




Western Blot analysis of POLR3F expression in transfected 293T cell line by POLR3F rabbit polyclonal antibody.

Lane 1: POLR3F transfected lysate(35.70 KDa).
Lane 2: Non-transfected lysate.

Antibody dilution: 1:1000. Western Blot analysis of POLR3F expression in HeLa (Predicted Mw: 35.6840 kDa).


Western Blotting (Cell lysate)



POLR3F rabbit polyclonal antibody. Western Blot analysis of POLR3F expression in HeLa.

Antibody dilution: 1:1000. Western Blot analysis of POLR3F expression in human liver (Predicted Mw: 35.6840 kDa).

Western Blotting (Tissue lysate)



POLR3F rabbit polyclonal antibody. Western Blot analysis of POLR3F expression in human liver.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619