

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.


Catalog No. CNP- WB-316

Anti SURF1

Product type	Primary antibodies
Immunogen	SULF1 (NP_055985.1, 1 a.a. ~ 871 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MKYSCCALVLAVLGTELLGSLCSTVRSRFRGRIQQERKNIRPNILVLTDDQDVELGSL QVMNKTRKIMEHGGATFINAFVTPMCCPSRSSMLTGKYVHNHNVTNNENCSSPSWQAM HEPRTFAVYLNNTGYRTAFFGKYLNEYNGSYIPPGWREWLGLIKNSRFYNYTVCNGIKE KHGFDYAKDYFTDLITNESINYFKMSKRMYPHRPVMVISHAAPHGPEDSAPQFSKLYPN ASQHITPSYNYAPNMDKHWIMQYTGPMPLIHMEFTNILQRKRLQTLMSVDDSDVERLYNML VETGELENTYIIYTADHGYHIGQFGLVKGKSMPYDFDIRVPFFIRGPSVEPGSIVPQIVL NIDLAPTILDIAGLDTPPDVDGKSVLKLLDPEKPGNRFRTNKKAKIWRDFTLVERGKFLR KKEESSKNIQQSNHLPKYERVKELCQQARYQTACEQPGQKWQCIEDTSGKLRIHKCKGPS DLLTVRQSTRNLYARGFHDKDKECSCRESGYRASRSQRKSQRQFLRNQGTPKYKPRFVHT RQTRSLSVFEGEYDINLEEEELQVLQPRNIAKRHDEGHKGPRDLQASSGGNRGRMLA DSSNAVGPPTTVRVTHKCFILPNDSIH CERELYQSARAWKDHKAYIDKEIEALQDKIKNL REVRGHLKRRKPEECSCSKQSYYNKEKGVKKQEKLKSHLHPFKEAAQEVD SKLQLFKENN RRRKKEKRRQRKGEECSLPGLTCFTHDNNHWQTAPFWNLGSFCACTSSNNNTYWCLR TVNETHNFLFCEFATGFLEYFDMNTDPYQLTNTVHTVERGILNQLHVQLMELRSCQGYKQ CNP RPKNLDVGNKDGGSYDLHRGQLWDGWEG
Gene Alias	FLJ30905, FLJ38022, FLJ41750, HSULF-1, KIAA1077, SULF-1

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~101.00 kDa using SULF1 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of SULF1 expression in transfected 293T cell line by SULF1 rabbit polyclonal antibody.

Lane 1: SULF1 transfected lysate(101.00 KDa).
Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619