

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-351

Anti CYP4F22

Product type	Primary antibodies
Immunogen	CYP4F22 (AAH69351.1, 1 a.a. ~ 531 a.a) full-length human protein
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MLPITDRLLHLLGLEKTA FRIYAVSTLLLFLFFLFRLLLRFLRLCRSFYITCRRRLRCFP QPPRRNWLLGHLGMYLPNEAGLQDEKKVLDNMHHVLLVWMGPVLP LLVLVHPDYIKPLL GASAAIAPKDDLFYGF LKPWLG DGLLSKGDKWSRHRLLTPAFHFDILKPYMKIFNQ SADIMHAKWRHLAEGSAVSLDMFEHISLMTLDSLQKCVFSYNSNCQEKMSDYISAI IELSALS VRRQYRLHHYLDFIYYRSADGRRFRQACDMVHHFTTEVIQERRRALRQQGAEAWLKAKQG KTLDFIDVLLLARDEDGKELSD EDIRAEADTFMFEGHDTTSSGISWMLFNLAKYPEYQEK CREEIQEVMKGRELEEELEWDDLTQLPFTTMCIKESLRQYPPVTLVSRQCTEDIKLPDGRI IPKGIICLVSIYGT HHNPTVWPDSKVYNPYRFDPDNPQQRSPLAYVVPFSAGPRNCIGQSF AMAELRVVALTLRFRLSVDRT RKVRRKPELILRTENGLWLKVEPLPPRA


FLJ39501, LI3

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~62.00 kDa using CYP4F22 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of CYP4F22 expression in transfected 293T cell line by CYP4F22 rabbit polyclonal antibody.

Lane 1: CYP4F22 transfected lysate(62.00 KDa).
Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619