

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-352

Anti CYP20A1

Product type Primary antibodies

Immunogen CYP20A1 (AAH20616.1, 1 a.a. ~ 462 a.a) full-length human protein.

Form Liquid
Protein A purified

Volume 100ul

Specificity Human

Host Rabbit

Immunogen Sequence MLDFAIFAVTFLLALVGAVLYLYPASRQAAGIPGITPTEEKDGNLPDIVNSGSLHEFLVN
LHERYGPVVSFWFGRRLVSLGTVDVLKQHINPNKTLDPFETMLKSLLRYQSGGGSVSEN
HMREKLYENGVTDSLKSNFALLKLSEELLDKWLSTYPETQHVPLSQHMLGFAMKSVTQMV
MGSTFEDDQEVIRFQKNHGTWVSEIGKGFLDGLDKNMTRKKQYEDALMQLESVLRNIK
ERKGRNFSQHIFIDSLVQGNLNDQQILED SMIFSLASCIITAKLCTWAICFLTTSSEEVQK
KLYEEINQVFGNGPVTPEKIEQLRYCQHVLCE TVRTAKLTPVSAQFQDIEGKIDRFIIPR
ETLVLYALGVVLQDPNTWPSPHKFD PDRFDDEL VMKTFSSLGFGSGTQCEPELRFAYMVT
VLLSVLVKRLHLLSVEGQVIETKYELVTSSREEAWITVSKRY


CYP-M, MGC22229

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~52.50 kDa using CYP20A1 transfected 293T lysate.

Western Blotting (Transfected lysate)




Western Blot analysis of CYP20A1 expression in transfected 293T cell line by CYP20A1 rabbit polyclonal antibody.

Lane 1: CYP20A1 transfected lysate(52.50 KDa).

Lane 2: Non-transfected lysate.

Antibody dilution: 1:1000. Western Blot analysis of CYP20A1 expression in HeLa (Predicted Mw: 52.4320 kDa).

Western Blotting (Cell lysate)



CYP20A1 rabbit polyclonal antibody. Western Blot analysis of CYP20A1 expression in HeLa.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619