



POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-359

Anti HABP4


Product type	Primary antibodies
Immunogen	HABP4 (AAH18788.1, 1 a.a. ~ 308 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MKGALGSPVAAAGAAMQESFGCVVANRFHQLLDDESDPFDILREAERRRQQQLQRKRRDE AAAAAGAGPRGGRSPAGASGHRAGAGGRRESQKERKSLPAPVAQRPDSPGGGLQAPGQKR TPRRGEQQGWNSDRGPEGMLERAERSYREYRPHYETERQADFTAEEKFPDEKVPPELEVEEE TQVQEMTLDEWKNLQEQTRPKPEFNIRKPESTVPSKAVVIHKSRYRDMVMKDDYEDDSHV FRKPANDITSQLEINFGNLPRPGRGARGGTRGGGRIRRAENYGPRAEVVMQDVAPNPDD PEDFPALS

IHABP4, Ki-1/57, MGC825, SERBP1L

Application notes

WB
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~34.40 kDa using HABP4 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of HABP4 expression in transfected 293T cell line by HABP4 rabbit polyclonal antibody.

Lane 1: HABP4 transfected lysate(34.40 KDa).
Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619