

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-368

Anti C12orf24

Product type	Primary antibodies
Immunogen	C12orf24 (AAH20967.1, 1 a.a. ~ 273 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MLGQLLPHTARGLGAAEMPGQGPGSDWTERSSSAEPPAVAGTEGGGGGSAGYSCYQNSKG SDRIKDGKYNVSHIAKLQELWKTPQNQTIHLSKSMMEASFFKHPDLTTGQKRYLCSIAKI YNANYLKMLMKRQYMHVLQHSSQKPGVLTHHRSRLSSRYSQKHYPCTTWRHQLEREDSG SSDIAAASAPEMLIQHSLWRPVRNKEGIKTGYASKTRCKSLKIFRRPRKLFMQTVSSDDS ESHMSEEKKEEDLLNNFMQSMSIEEQGEHMLT


HSU79274

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~30.80 kDa using C12orf24 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of C12orf24 expression in transfected 293T cell line by C12orf24 rabbit polyclonal antibody.

Lane 1: C12orf24 transfected lysate(30.80 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619