

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-386

Anti OSBPL10

Product type	Primary antibodies
Immunogen	OSBPL10 (NP_060254.2, 1 a.a. ~ 764 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit


Immunogen Sequence

MERAVQGTGDDGGGSSNSSSRSSSRATSAGSSPSCSLAGRGVSSRSAAAGLGGGGSRSSPGS
VAASPSGGGGRRREPALLEGVLSKYTNLLQGWQNRYFVLDFEAGILQYFVNEQSKHQKPRG
VLSLSGAIVSLSDDEAPHMLVVYSANGEMFKLRAADAKEKQFWVTQLRACAKYHMEMNSKS
APSSRSRSLTLLPHGTPNSASPCSQRHLSVGAPGVVTITHHKSPAAARRAKSQYSGQLHE
VREMMNQVEGQKQNLVHAIESLPGSGPLTALDQDLLLLKATSAATLSCLGECLNLLQQSV
HQAGQPSQKPGASENIGWHGSKSHSTEQLKNGTLGSLPSASANITWAILPNSAEDEQTS
QPEPEPNSGSELVLSDEKSDNEDKEETELGVMEDQRSIILHLISQLKLGMDLTKVVLPT
FILEKRSLLMEYADFMHPDLLLLAITAGATPEERVICFVEYYLTAFHEGRKGALAKKPYN
PIIGETFHCSEWVPKDRVKPKRTASRSPASCHEHPMADDPKSKYKLRFAEQVSHHPPIS
CFYCECEEKRLCVNTHVWTKSKFMGMSVGVSMIGEVLRLLHGEYVFTLPSAYARSIL
TIPWVELGGKVSINCAKTGYSATVIFHTKPFYGGKVHRVTAEVKHNPTNTIVCKAHGEWN
GTLEFTYNNGETKVIDTTTTLPVYPKKIRPLEKQGPMEARNLWREVTYRLRLGDIDAATEQ
KRHLEEKQRVEERKRENLRTPWPKYFIQEGDGWVYFNPLWKAH

FLJ20363, ORP10, OSBP9

Application notes Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~84.00 kDa using OSBPL10 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of OSBPL10 expression in transfected 293T cell line by OSBPL10 rabbit polyclonal antibody.

Lane 1: OSBPL10 transfected lysate(84.00 KDa).
Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619