



POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-389

Anti GINS1


Product type	Primary antibodies
Immunogen	GINS1 (NP_066545.2, 1 a.a. ~ 196 a.a) full-length human protein.
Form	Liquid
	Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit

Immunogen Sequence MFCEKAMELIRELHRAPEGQLPAFNEDGLRQVLEEMKALYEQNQSDVNEAKSGGRSDLIP
TIKFRHCSSLRNRRCTVAYLYDRLLRIRALRWEYGSILPNALRFHMAAEEMEWFNYYKRS
LATYMRSLLGGDEGLDITQDMKPPKSLYIEVRCLKDYGEFEVDDGTSVLLKKNSQHFLPRW
KCEQLIRQGVLEHLS

Gene Alias KIAA0186, PSF1

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~23.00 kDa using GINS1 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of GINS1 expression in transfected 293T cell line by GINS1 rabbit polyclonal antibody.


Lane 1: GINS1 transfected lysate(23.00 kDa).

Lane 2: Non-transfected lysate.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of GINS1 transfected lysate.

Immunoprecipitation



Immunoprecipitation of GINS1 transfected lysate using GINS1 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with GINS1 purified mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619