

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-395

Anti SAMHD1

Product type	Primary antibodies
Immunogen	DCPS (NP_054745.1, 1 a.a. ~ 337 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit

Immunogen Sequence

MADAAPQLGKRKRELDVEEAHAASTEKEAGVGNGTCAVRLPFSGFRLQKVLRESARDK
IIFLHGKVNEASGDGDGEDAVVILEKTPFQVEQVAQLLTGSPQLQFSNDIYSTYHLFP
PRQLNDVKTTVVYPATEKHLQKYLQDLRLIRETGDDYRNITLPHLESQSLSIQWVYNIL
DKKAEADRIVFENPDPSDGFVLIPDLKWNQQQLDDLYIAICHRRGIRSLRDLTPEHLPL
LRNILHQQQEAILQRYRMKGDHLRVYLHYLPSYYHLHVHFTALGFEAPGSGVERAHLAE
VIENLECDPRHYQQRTLTFALRADDPLLKLLQEAQQS


HINT-5, HSL1, HSPC015

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~38.60 kDa using DCPS transfected 293T lysate

Western Blotting (Transfected lysate)




Western Blot analysis of DCPS expression in transfected 293T cell line by DCPS rabbit polyclonal antibody.

Lane 1: DCPS transfected lysate(38.60 KDa).

Lane 2: Non-transfected lysate.

Antibody dilution: 1:1000. Western Blot analysis of DCPS expression in mouse liver (Predicted Mw: 38.6090 kDa).

Western Blotting (Tissue lysate)



DCPS rabbit polyclonal antibody. Western Blot analysis of DCPS expression in mouse liver.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619