

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-400

Anti CTPS2

Product type Primary antibodies

Immunogen CTPS2 (NP_062831.3, 1 a.a. ~ 586 a.a) full-length human protein.

Form Liquid
Protein A purified

Volume 100ul

Specificity Human

Host Rabbit

Immunogen Sequence MKYILVTGGVISGIGKGIASSIGTILKSCGLRVTAIKIDPYINIDAGTFSPYEHGEV
FV LNDGGEVDLDLGNRYERFLDINLYKDNNITGKIYQHVINKERRGDYLGKTVQVPHITDA
VQEWVMNQAKVPVDGNKEEPQICVIELGGTIGDIEGMPFVEAFRQFQFKAKRENFCNIHV
SLVPQLSATGEQKTKPTQNSVRALRGLGLSPDLIVCRSSTPIEMAVKEKISMFCNVNPEQ
VICIHDVSSTYRVPVLLLEEQSIVKYFKERLHLPIGDSASNLLFKWRNMADRYERLQKICS
IALVGKYTKLRDCYASVFKALEHSALAINHKLNLMYIDSIDLEKITETEDPVKFHEAWQK
LCKADGILVPGGFGIRGTLGKLQAISWARTKKIPFLGVCLGMQLAVIEFARNCLNLKDDAD
STEFRPNAPVPLVIDMPEHNPGNLGGTMRLGIRRTVFKTENSILRKLYGDVPFIEERHRH
RFEVNPNIKQFEQNDLSFVGQDVGDRMEIIELANHPYFVGVQFHPEFSSRPMKPSPPY
LGLLLAATGNLNAYLQQGCKLSSSDRYSDASDDSFSEPRIAELEIS


DKFZp686C17207, FLJ43358, MGC32997

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~65.70 kDa using CTPS2 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of CTPS2 expression in transfected 293T cell line by CTPS2 rabbit polyclonal antibody.

Lane 1: CTPS2 transfected lysate(65.70 kDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619