

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.


Catalog No. CNP- WB-404

Anti FAM117A

Product type	Primary antibodies
Immunogen	FAM117A (NP_110429.1, 1 a.a. ~ 453 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MAGAAAGGRGGGAWGPGRGGAGGLRRGCSPAPAGSPRAGLQPLRATIPFQLQQPHQRRD GGGRAASVPCSVAPKSVCRPQPLQVRRTFSLDTILSSYLLGQWPRDADGAFTCCTNDKA TQTPLSWQELEGERRASSCAHKRSASWGSTDHRKEISKLKQQLQRTKLSRSGKEKERGSPL LGDHAVRGALRASPPSFPSGSPVLRRLSPCLHRSLEGLNQEEVFVKEQGEEELLRLDI PDGHRAPAPPQSGSCDHPLLLLEPGNLASSPSMSLASPQPCGLASHEEHARGAAEELASTP NDKASSPGHPAFLEDGSPSPVLAFASPRPNHSYIFKREPPEGCEKVRVFEEATSPGPD AFLTSCPDKNKVHFNPTGSAFCPVNLMKPLFPGMGFIFRNCPSNPGSPLPPASPRPPRK DPEASKASPLPFEPWQRTTPPSEEPVLFQSSLMV

Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~48.30 kDa using FAM117A transfected 293T lysate.
--------------------------	---

Western Blotting (Transfected lysate)



Western Blot analysis of FAM117A expression in transfected 293T cell line by FAM117A rabbit polyclonal antibody.

Lane 1: FAM117A transfected lysate(48.30 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619