

**POLYCLONAL ANTIBODY**

*For research use only. Not for clinical diagnosis.*

**Catalog No. CNP- WB-407**

# Anti CXADR

**Product** Primary antibodies

**type**

**Immunogen** CXADR (NP\_001329.1, 1 a.a. ~ 365 a.a) full-length human protein.

**Form**

Liquid

Protein A purified

**Volume**

100ul

**Specificity**

Human

**Host**

Rabbit

**Immunogen**

**Sequence**

MALLLCFVLLCGVVDFAFSLSSITTPPEEMIEKAKGETAYLPCKFTLSPEDQGGLDIEWLIS  
PADNQKVDQVILYSGDKIYDDYYPDLKGRVHFTSNDLKSGDASINVTNLQLSDIGTYQC  
KVKKAPGVANKKIHLVVLVKPSGARCVDGSEEIGSDFKIKCEPKEGSLPLQYEWQKLSD  
SQKMPTSWLAEMTSSVISVKNASSEYSGTYSCTVRNRVGSQCLLRNLNVPPSNKAGLIA  
GAIIGTLLALALIGLIIFCCRKKRREEKYEKEVHHDIREDVPPPKSRTSTARSYIGSNHS  
SLGSMSPSNMEGYSKTQYNQVPSEDFERTPQSPTLPPAKVAAPNLSRMGAIPVMIPAQSK  
DGSIV

CAR, HCAR


**Application**

**notes**

**WB**

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~40.00 kDa using CXADR transfected 293T lysate.

**Western Blotting (Transfected lysate)**




Western Blot analysis of CXADR expression in transfected 293T cell line by CXADR rabbit polyclonal antibody.

Lane 1: CXADR transfected lysate(40.00 KDa).

Lane 2: Non-transfected lysate.

Antibody dilution: 1:1000. Western Blot analysis of CXADR expression in Jurkat (Predicted Mw: 40.0300 kDa).

**Western Blotting (Cell lysate)**



CXADR rabbit polyclonal antibody. Western Blot analysis of CXADR expression in Jurkat.

**Storage**

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

*For research use only. Not for clinical diagnosis.*



**COSMO BIO Co., LTD.**

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp)

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619