

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-472

Anti CDC34

Product type Primary antibodies

Immunogen CDC34 (NP_004350.1, 1 a.a. ~ 236 a.a) full-length human protein.

Form Liquid
Protein A purified

Volume 100ul

Specificity Human


Host Rabbit

Immunogen Sequence MARPLVPSSQKALLLELKGLEQEEPVGEFRVTLVDEGDLYNWEVAIFGPPNTYYEGGYFKA
RLKFPIDYPYSPPAFRFLTKMWHPNYETGDVCISILHPPVDDPQSGELPSERWNPTQNV
RTILLSVISLLNEPNTFSPANVDASVMYRKWKESKGDREYTDIIRKQVLGTVDAERDG
VKVPTTLAEYCVKTKAPAPDEGSDLFYDDYYEDGEVEEEEADSCFGDDEDDSGTEES

Gene Alias E2-CDC34, UBC3, UBE2R1

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~26.70 kDa using CDC34 transfected 293T lysate..

Western Blotting (Transfected lysate)



Western Blot analysis of CDC34 expression in transfected 293T cell line by CDC34 rabbit polyclonal antibody.

Lane 1: CDC34 transfected lysate(26.70 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619