

POLYCLONAL ANTIBODY


For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-489

Anti GPS2

Product type	Primary antibodies
Immunogen	GPS2 (NP_004480.1, 1 a.a. ~ 327 a.a) full-length human protein
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MPALLERPKLSNAMARALHRHIMMERERKRQEEEEVDKMMEQKMKEEQERRKKKEMEERM SLEETKEQILKLEEKLLALQEEKHQFLQLKKVLHEEEKRRRKEQSDLTTLTSAAYQQSL TVHTGTHLLSMQGSPPGHNRPGLMAADRAKQMFQPQVLTTTRHYVGSAAAFAGTPEHGQF QGSPGGAYGTAQPPPHYGPTQPAYSPSQQLRAPSAFPAVQYLSQPQPQPYAVHGHFQPTQ TGFLQPGGALSLLQKQMEHANQQTGFSDSSSLRPMHPQALHPAPGLLASPQLPVQMMPAGK SGFAATSQPGPRLPFIQHSQNPRFYHK
Gene Alias	AMF-1, MGC104294, MGC119287, MGC119288, MGC119289
Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~36.70 kDa using GPS2 transfected 293T lysate

Western Blotting (Transfected lysate)



Western Blot analysis of GPS2 expression in transfected 293T cell line by GPS2 rabbit polyclonal antibody.

Lane 1: GPS2 transfected lysate(36.70 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619