

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-507

Anti RBBP4

Product type Primary antibodies

Immunogen RBBP4 (NP_005601.1, 1 a.a. ~ 425 a.a) full-length human protein

Form Liquid
Protein A purified

Volume 100ul

Specificity Human

Host Rabbit

Immunogen Sequence MADKEAAFDDAVEERVINEEYKIWKKNTPFLYDLVMTHALEWPSLTAQWLDPDTRPEGKD
FSIHLVLGTHTSDEQNHLVIASVQLPNDDAQFDASHYDSEKGEFGGFGSVSGKIEIEIK
INHEGEVNRARYMPQNPCIATKTPSSDVLVFDYTKHPSKPDPSGECNPDRLRLRGHQKEG
YGLSWNPNLSGHLLSASDDHTICLWDISAVPKEGKVVDKTIPTGHTAVVEDVSWHLLHE
SLFGSVADDQKLMIWDTRSNNSTKPSHSDAHTAEVNCLSFNPYSEFILATGSADKTVAL
WDLRNLKLLHFSFESHKDEIFQVQWSPHNETILASSGTDRRLNVWDLKIGEEQSPEDAE
DGPPELLFIHGGHTAKISDFS WNPNEPWVICSVSEDNIMQVWQMAENIYNDEDPEGSVDP
EGQGS


Gene Alias NURF55, RBAP48

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~47.70 kDa using RBBP4 transfected 293T lysate.

Western Blotting (Transfected lysate)




Western Blot analysis of RBBP4 expression in transfected 293T cell line by RBBP4 rabbit polyclonal antibody.

Lane 1: RBBP4 transfected lysate(47.70 KDa).

Lane 2: Non-transfected lysate.

Antibody dilution: 1:1000. Western Blot analysis of RBBP4 expression in HeLa (Predicted Mw: 47.6560 kDa).

Western Blotting (Cell lysate)



RBBP4 rabbit polyclonal antibody. Western Blot analysis of RBBP4 expression in HeLa.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619