

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-517

Anti TP53

Product type Primary antibodies

Immunogen TP53 (NP_000537.2, 1 a.a. ~ 393 a.a) full-length human protein.

Form Liquid
Protein A purified

Volume 100ul

Specificity Human

Host Rabbit

Immunogen Sequence MEEPQSDPSVEPPLSQETFSDLWKLLPENNVLSPLPSQAMDDLMLSPDDIEQWFTEDPGP
DEAPRMPEAAPRVAPAPAAPTPAAPAPAPSWPLSSSVPSQKTYQGSYGFRGLGFLHSGTAK
SVTCTYSPALNKMFCQLAKTCPVQLWVDSTPPPGRVVRAMAIYKQSQHMTEVVRRCPPHHE
RCSDSDGLAPPQHILIRVEGNLRVEYLDDRNTFRHSVVVPYEPPEVGSDCTTIHYNMCNS
SCMGGMNRRPILTIITLEDSSGNLLGRNSFEVRVCACPGRRRTEENLRKKGEPHHELP
PGSTKRALPNNTSSSPQPKKKPLDGEYFTLQIRGRERFEMFRELNEALELKDAQAGKEPG
GSRAHSSHLKSKKGQSTSRHKKLMFKTEGPDSD


Gene Alias FLJ92943, LFS1, TRP53, p53

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~43.23 kDa using TP53 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of TP53 expression in transfected 293T cell line by TP53 rabbit polyclonal antibody.

Lane 1: TP53 transfected lysate(43.23 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619