

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-526

Anti ACOX3

Product type Primary antibodies

Immunogen ACOX3 (AAH17053.1, 1 a.a. ~ 624 a.a) full-length human protein.

Form Liquid
Protein A purified

Volume 100ul

Specificity Human

Host Rabbit

Immunogen Sequence MASTVEGGDTALLPEFPRGPLDAYRARASFSWKELALFTEGEGMLRFKKTIFSALENDPL
FARSPGADLSLEKYRELNFLRCKRIFEYDFLSVEDMFKSPLKVPALIQCLGMYDSSLA
AKYLLHSLVFGSAVYSSGSRHLYTIQKIFRMEIFGCFALTELSHGSNTKAIRTTAHYDPAT
EEFIHSPDFEAAKFWVGNMGKTATHAVVFAKLCVPGDQCHGLHPFIVQIRDPKTLLPMP
GVMVGDIGKKLGQNGLDNGFAMFHKVRVPRQSLNRMGDVTPEGTYVSPFKDVRQRF
GASLGSLSSGRVSIVSLAILNLKLAVALRFSATRRQFGPTEEEEIPVLEYPMQQRLLPYL
AAVYALDHFSKSLFLDLVELQRGLASGDRSARQAEIGREIHALASASKPLASWTTQGG
IQECREACGGHGYLAMNRLGVLRDDNDPNCTYEGDNNILLQQTSNYLLGLLAHQVHDG
ACFRSPLKSVDFLDAYPGILDQKFEVSSVADCLDSAVALAAYKWLVCYLLRETYQKL
NQEKRSGLSSDFEARNKCQVSHGRPLALAFVELTVVQRFHEHVVHQPSPVPSLRAVL
GRLSALYALWSLSRHAALLYRAERRCSCPGRDRSS


Gene Alias

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~69.60 kDa using ACOX3 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of ACOX3 expression in transfected 293T cell line by ACOX3 rabbit polyclonal antibody.

Lane 1: ACOX3 transfected lysate(69.60 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619