

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-528

Anti DDO

Product type Primary antibodies

Immunogen DDO (NP_003640.2, 1 a.a. ~ 369 a.a) full-length human protein.

Form Liquid
Protein A purified

Volume 100ul

Specificity Human

Host Rabbit

Immunogen Sequence MRPARHWETRFGARDFGGFQDCFFRDRLMDTARIAVVGAGVVGLSTAVCISKLVPRCSVT
IISDKFTPDTTSDVAAGMLIPHTYPDTPIHTQKQWFRET FNHLFAIANSAEAGDAGVHLV
SGWQIFQSTPTEEVFWADVVLGFRKMTEAELKKFPQYVFGQAFTTLKCECPAYLPWLEK
RIKSGSGWTLTRRIEDLWELHPSFDIVVNCSGLGSRQLAGDSKIFPVRGQVLQVQAPWVE
HFIRDGSLTYIYPGTSHVTLGGTRQKGDWNLSPDAENSREILSRCCALEPSLHGACNIR
EKVGLRPYRPGVRLQTELLARDGQRLPVVHHYGHGSGGISVHWGTALEAARLVSECVHAL
RTPIPKSNL


Gene Alias DASOX, DDO-1, DDO-2, FLJ45203

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~41.00 kDa using DDO transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of DDO expression in transfected 293T cell line by DDO rabbit polyclonal antibody.

Lane 1: DDO transfected lysate(41.00 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619