



POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-530

Anti FADS2

Product type Primary antibodies

Immunogen FADS2 (NP_004256.1, 1 a.a. ~ 444 a.a) full-length human protein

Form Liquid
Protein A purified

Volume 100ul

Specificity Human

Host Rabbit

Immunogen Sequence MGKGGNQGEAAEREVSPTFSWEEIQKHNLRTDRWLVIDRKVYNITKWSIQHPGGQQRVI
GHYAGEDATDAFRAHPDLEFVGKFLKPLLIGELAPEEPSQDHGKNSKITEDFRALRKTA
EDMNLFKTNHVFFLLLAHIIALESIAWFTVFYFGNGWIPTLITAFVLATSQAQAGWLQH
DYGHLSVYRKPKWNHLVHKFVIGHLKGASANWWNHRHFQHHAKPNIFHKDPDVNMLHVFV
LGEWQPIEYGKKKLYLPYNHQHEYFFLIGPPLLIPMYFQYQIIMTMIVHKNWVDLAWAV
SYYIRFFITYIPFYGILGALLFLNFIRFLESHWVWVTQMNHIVMEIDQEAYRDWFSSQL
TATCNVEQSFFNDWFSGHLNFQIEHHLFPTMPRHNHLHKIAPLVKSLCAKHGIEYQEKPLL
RALLDIIRSLKSGKLWLDAYLHK


Gene Alias D6D, DES6, FADSD6, LLCDL2, SLL0262, TU13

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~52.30 kDa using FADS2 transfected 293T lysate..

Western Blotting (Transfected lysate)



Western Blot analysis of FADS2 expression in transfected 293T cell line by FADS2 rabbit polyclonal antibody.

Lane 1: FADS2 transfected lysate(52.30 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619