

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- WB-533

Anti ARM CX2

Product type Primary antibodies

Immunogen ARM CX2 (NP_055597.1, 1 a.a. ~ 632 a.a) full-length human protein

Form Liquid
Protein A purified

Volume 100ul

Specificity Human

Host Rabbit

Immunogen Sequence MSRVRDAGCVAAGIVIGAGAWYCVYKYTRGRDQTKKRM AKPKNRAVAGTGARARAGLRAG
FTIDLGS GFSPPTPVRAEAEDRAQDEASALDTVGA EAVAPAASSAE AQSGAGSQAQ EADG
AGVGPKAESV VGAAMASAIAPPPGVTEALGAAEAPAMAGAPKVAEAPREAETSRAAVPPG
TVVPTEAAAPTEVTEGPGVAAPTKVAEAPGVASPT EAAEAPVPATPTGAAAPTGA AESPG
TSGSPRTAVVPGTSAKKATPGAHTGAIPKATSATGAVPKGGGKG VTRSRNGGKGKGKKS
KVEVDELGMGFRPGDGA AAAAAAASANGGQAFLAEVDPDSEEGESGWTDTESDSDSEPETQR
RGRGRRPVAMQKRPFPEIDEILGVRDLRKVLALLQKSDDPFIQQVALLTLSNNANYSCN
QETIRKLGLPIIANMINKTDPHIKEKALMAMNNLSENYENQGR LQVYMKNKVMDDIMASN
LNSAVQVVGLKFLT NMTITNDYQHLLVNSIANFFRLLSQGGGKIKVEILKILSNFAENPD
MLKKLLSTQVPASFSSLYNSYVESEILINALTFEIIYDNLRAEVFN YREFNKGSLFYLC
TTSGVCVKKIRALANHHDLLVKVKVIKLVNKF


Gene Alias ALEX2, KIAA0512, MGC13343, MGC8742

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~65.70 kDa using ARM CX2 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of ARM CX2 expression in transfected 293T cell line by ARM CX2 rabbit polyclonal antibody.

Lane 1: ARM CX2 transfected lysate(65.70 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619